

IME 1050.5 SERIES

Monolithic un-packaged
 Monolithic-packaged
 Combined Differential-un-packaged
 Combined Differential-packaged
 Full Differential-un-packaged
 Full Differential-packaged

IME 0550.5 SERIES

Monolithic un-packaged
 Monolithic-packaged
 Combined Differential-un-packaged
 Combined Differential-packaged
 Full Differential-un-packaged
 Full Differential-packaged

IME 2025.3 SERIES

Monolithic un-packaged
 Monolithic-packaged
 Full Differential-un-packaged
 Full Differential-packaged

IME 1525.3 SERIES

Monolithic un-packaged
 Monolithic-packaged
 Full Differential-un-packaged
 Full Differential-packaged

IME 1025.3 SERIES

Monolithic un-packaged
 Monolithic-packaged
 Full Differential-un-packaged
 Full Differential-packaged

IME 0525.3 SERIES

Monolithic un-packaged
 Monolithic-packaged
 Full Differential-un-packaged
 Full Differential-packaged

GOLD-Au

IME 1050.5-M-Au
 IME 1050.5-M-Au-P
 IME 1050.5-CD-Au-U
 IME 1050.5-CD-Au-P
 IME 1050.5-FD-Au-U
 IME 1050.5-FD-Au-P

GOLD-Au

IME 0550.5-M-Au
 IME 0550.5-M-Au-P
 IME 0550.5-CD-Au-U
 IME 0550.5-CD-Au-P
 IME 0550.5-FD-Au-U
 IME 0550.5-FD-Au-P

GOLD-Au

IME 2025.3-M-Au
 IME 2025.3-M-Au-P
 IME 2025.3-FD-Au-U
 IME 2025.3-FD-Au-P

GOLD-Au

IME 1525.3-M-Au
 IME 1525.3-M-Au-P
 IME 1525.3-FD-Au-U
 IME 1525.3-FD-Au-P

GOLD-Au

IME 1025.3-M-Au
 IME 1025.3-M-Au-P
 IME 1025.3-FD-Au-U
 IME 1025.3-FD-Au-P

GOLD-Au

IME 0525.3-M-Au
 IME 0525.3-M-Au-P
 IME 0525.3-FD-Au-U
 IME 0525.3-FD-Au-P

PLATINUM-Pt

IME 1050.5-M-Pt-U
 IME 1050.5-M-Pt-P
 IME 1050.5-CD-Pt-U
 IME 1050.5-CD-Pt-P
 IME 1050.5-FD-Pt-U
 IME 1050.5-FD-Pt-P

PLATINUM-Pt

IME 0550.5-M-Pt-U
 IME 0550.5-M-Pt-P
 IME 0550.5-CD-Pt-U
 IME 0550.5-CD-Pt-P
 IME 0550.5-FD-Pt-U
 IME 0550.5-FD-Pt-P

PLATINUM-Pt

IME 2025.3-M-Pt-U
 IME 2025.3-M-Pt-P
 IME 2025.3-FD-Pt-U
 IME 2025.3-FD-Pt-P

PLATINUM-Pt

IME 1525.3-M-Pt-U
 IME 1525.3-M-Pt-P
 IME 1525.3-FD-Pt-U
 IME 1525.3-FD-Pt-P

PLATINUM-Pt

IME 1025.3-M-Pt-U
 IME 1025.3-M-Pt-P
 IME 1025.3-FD-Pt-U
 IME 1025.3-FD-Pt-P

PLATINUM-Pt

IME 0525.3-M-Pt-U
 IME 0525.3-M-Pt-P
 IME 0525.3-FD-Pt-U
 IME 0525.3-FD-Pt-P

INDIUM TIN OXIDE-ITO

IME 1050.5-M-ITO-U
 IME 1050.5-CD-ITO-U
 IME 1050.5-FD-ITO

INDIUM TIN OXIDE-ITO

IME 0550.5-M-ITO-U
 IME 0550.5-CD-ITO-U
 IME 0550.5-FD-ITO

INDIUM TIN OXIDE-ITO

IME 2025.3-M-ITO-U
 IME 2025.3-FD-ITO

INDIUM TIN OXIDE-ITO

IME 1525.3-M-ITO-U
 IME 1525.3-FD-ITO

INDIUM TIN OXIDE-ITO

IME 1025.3-M-ITO-U
 IME 1025.3-FD-ITO

INDIUM TIN OXIDE-ITO

IME 0525.3-M-ITO-U
 IME 0525.3-FD-ITO

INTERDIGITATED MICROSENSOR ELECTRODES (IMES)

| IME DEVICES | Line and Space, Digit Length | Designs | Conductor |
|--------------------------|------------------------------|------------------|--------------------|
| IME 2050.5 SERIES | 20 microns, 5 mm long | M | Au, Pt, ITO |
| IME 1550.5 SERIES | 15 microns, 5 mm long | M, CD, FD | Au, Pt, ITO |
| IME 1050.5 SERIES | 10 microns, 5 mm long | M, CD, FD | Au, Pt, ITO |
| IME 0550.5 SERIES | 5 microns, 5 mm long | M, CD, FD | Au, Pt, ITO |
| IME 2025.3 SERIES | 20 microns, 3 mm long | M, CD, FD | Au, Pt, ITO |
| IME 1525.3 SERIES | 15 microns, 3 mm long | M, CD, FD | Au, Pt, ITO |
| IME 1025.3 SERIES | 10 microns, 3 mm long | M, CD, FD | Au, Pt, ITO |
| IME 0525.3 SERIES | 5 microns, 3 mm long | M, CD, FD | Au, Pt, ITO |

ABTECH -- Chemical and biological sensor devices, instruments, and sensor systems.

■ **Interdigitated Microsensor Electrodes (IMEs)** are inert, array microelectrodes formed from microlithographically patterned conductors onto an insulating substrate chip. They are designed for the simultaneous interrogation of the electrical, electrochemical, and optical properties of thin polymeric films and coatings, for applications in microelectrochemistry, for electrical/electrochemical impedance spectroscopy, and for chemical and biological sensor development. ■ **Microfabricated** from magnetron sputter-deposited gold (Au), e-gun vapor-deposited platinum (Pt), or indium tin oxide (ITO) these devices occur in three configurations; Monolith (M), Combined Differential (CD) and Full Differential (FD), and as packaged electrodes (with attached leadwires and encapsulated) or as un-packaged chips. ■ **IME chips** are available with 5 μm, 10 μm, 15 μm or 20 μm line and space dimensions and of defined cell constant. ■ **Investigate** the chemoresistive responses of transducer-active, polymeric films in the same electrode configuration, the same test environment, and on the same sample film. ■ **In research and product development**, these devices are widely used for conductimetric/impedimetric, chemoresistive chemical and biological sensors using electrically conducting (electroconductive) polymers, for impedance sensors based on Langmuir-Blodgett thin films, for studying the environmental effects on polymer thin films, in micro-electrochemistry and in cell-based biochips. ■ **Develop** these devices into products where the application requires a compact, durable and versatile chemical or biological chemoresistive sensor of low cost.

At the forefront of Molecular Bioelectronics



For Further Information
 Customer Service and Support
 Tel.: + 1804.783.7829
 Fax: + 1804.783.7830
 URL: <http://www.abtechsci.com>
 Email: abtech@abtechsci.com

Advanced Biosensor and Biochip Technology 800 East Leigh Street, Richmond, Virginia 23219, U.S.A.

■ Introduction

Interdigitated Microsensor Electrodes (IMEs) is the registered trade name for a family of devices developed by ABTECH. Also called interdigitated arrays (IDAs), these devices are microfabricated (using microelectronics fabrication techniques) from patterns of noble metals deposited on an insulating substrate chip. They are designed for the simultaneous interrogation of the electrical, electrochemical

and optical properties of thin polymeric films and coatings, for applications in microelectrochemistry, and for electrical/electrochemical impedance spectroscopy.

■ Applications

Applications of IMEs in research and product development include:
 3 Conductimetric chemical sensors and biosensors¹ based on electroactive

polymers². These devices exploit the very large change in electrical impedance that accompanies oxidation/reduction of these polymer films.

3 Electrical and Electrochemical Impedance Spectroscopy of organic thin films and coatings³

3 Capacitance probes and humidity sensors, e.g. based on Langmuir-Blodgett films⁴.

3 Modern microelectrochemistry. The performance of electroanalysis in high impedance environments⁵.

3 Conductimetry. To determine the conductance of low conductivity media and conductimetric titrations.

■ Coatings

Film or coating application to the IME device may be achieved by dip coating, spin casting, spray painting, air-brushing, brush painting, by Langmuir-Blodgett thin film deposition, by electropolymerization, and by molecular self assembly

For further information, request Application Notes:

IME1 - *Interdigitated Microsensor Electrodes: Applications and References.*

IME2 - *Conductimetric Urea Biosensor Formed From Interdigitated Microsensor Electrodes*

■ Technical Specifications

| Substrate: | Schott D263 Borosilicate Glass | | | |
|--|--------------------------------|-------------------------------------|---------------------|---|
| Dielectric Constant, ϵ_r , at 1 MHz | | 6.7 | | |
| Dielectric Loss Angle, $\tan \delta$, at 1 MHz | | 61×10^{-4} | | |
| Electrical Resistivity (50 Hz) (250 °C) | | $1.6 \times 10^8 \Omega \text{ cm}$ | | |
| Coefficient of Linear Thermal Expansion α , 20-300 °C | | $7.2 \times 10^{-6} \text{ K}^{-1}$ | | |
| Refractive Index at 20°C, n_e ($\lambda = 546.1 \text{ nm}$) | | 1.5249 | | |
| Metallization: | 100 Å Ti /W 1,000 Å Au or Pt | | | |
| XX50 Series | | 2050.5 | 1550.5 | 1050.5 0550.5 |
| Digit length, d, (μm): | | 4,980 μm | 4,985 μm | 4,990 μm 4,995 μm |
| No. of digit pairs, N | | 50 | 50 | 50 50 |
| Digit Width, a, (μm): | | 20 μm | 15 μm | 10 μm 05 μm |
| Interdigit Space, a, (μm): | | 20 μm | 15 μm | 10 μm 05 μm |
| Spatial Periodicity, λ , (μm) | | 80 μm | 60 μm | 40 μm 20 μm |
| Zaretsky ^{6,7} Meander Length, M, (cm) | | 24.90 | 24.93 | 24.95 24.98 |
| Center Line or Serpentine Length ⁹ (cm) | | 49.70 | 49.65 | 49.60 49.55 |
| Cell Constant ⁸ (cm^{-1}) | | 0.040 | 0.040 | 0.040 0.040 |
| XX25 Series | | 2025.3 | 1525.3 | 1025.3 0525.3 |
| Digit length, d, (μm): | | 2,980 μm | 2,985 μm | 2,990 μm 2,995 μm |
| No. of digit pairs, N | | 25 | 25 | 25 25 |
| Digit Width, a, (μm): | | 20 μm | 15 μm | 10 μm 05 μm |
| Interdigit Space, a, (μm): | | 20 μm | 15 μm | 10 μm 05 μm |
| Spatial Periodicity, λ , (μm) | | 80 μm | 60 μm | 40 μm 20 μm |
| Zaretsky ^{6,7} Meander Length, M, (cm) | | 7.45 | 7.46 | 7.48 7.49 |
| Center Line or Serpentine Length (cm) | | 14.80 | 14.77 | 14.75 14.73 |
| Cell Constant ⁸ (cm^{-1}) | | 0.040 | 0.040 | 0.040 0.040 |


| IME Chip Dimensions | Un-packaged Die ($l \times w \times t$) | Packaged Electrode* ($l \times w \times t$) |
|---------------------|---|---|
| M, CD and FD | 2.00 x 1.00 x 0.05 cm | 13.2 x 1.38 x 0.7 cm |
| XX25.3 M | 1.00 x 0.50 x 0.05 cm | 13.2 x 1.38 x 0.7 cm |

- *Electrode Body: PVC-jacketed printed circuit board
- *Encapsulant: Epoxy header. Polyimide packaged chip.
- *Leadwires: Color coded, 30AWG stranded copper, shielded, and PVC jacketed.

■ References and Notes

1. (a) Lawrence, A "Conductimetric Enzyme Assays." *European Journal of Biochemistry* **1981**, *18*, pp. 221-225. (b) Lawrence A.; Moores, G. R. "Conductimetry in Enzyme Studies." *Europ. J. of Biochem.* **1972**, *24*, pp. 538 - 546. (c) Cullen, D. C.; Sethi, R. S.; Lowe, C. R. "Multi-Analyte Miniature Conductance Biosensor," *Analytica Chimica Acta* **1990**, *231*, pp. 33 - 40. (d) Mikkelsen, S. R.; Rechnitz, G.A. "Conductometric Transducers for Enzyme-Based Biosensors," *Anal. Chem.* **1989**, *61*(15), pp. 1737 - 1741. (e) Hanss, M.; Rey, A. "Application de la Conductimétrie a l'étude des Réactions Enzymatiques, Systeme Urée-Uréase," *Biochim. Biophys. Acta.* **1971**, *277*, pp. 630 - 638. (f) Norman F. Sheppard, Jr. and Anthony Guiseppi-Elie "Enzyme Sensors Based on Conductimetric Measurement"; *In Protocols in Biosensors Research*; Ashok Mulchandani and Kim R. Rogers, Editors; Humana Press, Totowa, NJ. (in press). (g) Anthony Guiseppi-Elie, Gordon G. Wallace, and Tomakazu Matsue "Chemical and Biological Sensors Based on Electrically Conducting Polymers" *In Handbook of Conductive Polymers* 2nd Edition (1998), T. Skotheim, R. Elsenbaumer and J. R. Reynolds Eds; Marcel Dekker, New York. (h) Norman F. Sheppard, Jr., David J. Mears, and Anthony Guiseppi-Elie "Model of a Conductimetric Urea Biosensor" *Biosensors and Bioelectronics* (1996), Vol. 11(10) 967 - 979.
2. (a) Anthony Guiseppi-Elie, James M. Tour, David L. Allara and Norman F. Sheppard, Jr. "Bioactive Polypyrrole Thin Films with Conductimetric Response to Analyte" *In, Electrical, Optical, and Magnetic Properties of Organic Solid State Materials*, Eds. A. K.-Y. Jen, C. Y.-C. Lee, L. R.alton, M. F. Rubner, G. E. Wnek, L. Y. Chiang, Mat. Res. Soc. Symp. Proc. Vol. 413; Materials Research Society, Pittsburgh, **1996**, p 439- 444. (b) A. Guiseppi-Elie, A. M. Wilson, J. M. Tour, T. W. Brockmann, P. Zhang, D. L. Allara "Specific Immobilization of Electropolymerized Polypyrrole Thin Films onto Interdigitated Microsensor Array Electrodes" *Langmuir* (1995), *11*(45), 1768.
3. Naoi, K; Ueyama, K.; Osaka, T. "Impedance Analysis of Ionic Transport in Polypyrrole-Polyazulene Copolymer and its Charge-Discharge Characteristics" *J. Electrochem. Soc.* **1990**, *137*(2) 494.
4. Kowel, S. T.; Zhou, G.-G.; Srinivasan, M. P.; Stroeve, P.; Higgins, B. G. "On-line Diagnostics for Langmuir-Blodgett Film Growth" *Thin Solid Films* **1985**, *134*, 209.
5. Bard, A. J.; Crayston, J. A.; Kittlesen, G. P.; Shea, T. V.; Wrighton, M. S. "Digital Simulation of the measured Electrochemical Response of Reversible Redox Couples at Microelectrode Arrays: Consequences Arising from Closely Spaced Ultramicroelectrodes" *Anal. Chem.* **1986**, *58*(11), 2321.
6. Zaretsky, M. C.; Mouyad, L.; Melcher, J. R. *IEEE Trans. Electr. Insul.* **1988**, *23*, 897.
7. The Zaretsky convention defines the meander length; $M = N \cdot d$
8. Sheppard, N. F.; Tucker, R. C.; Wu, C. "Electrical Conductivity Measurements Using Microfabricated Interdigitated Electrodes" *Anal. Chem.* **1993**, *65*, 1199.
9. Serpentine length is defined as: $S = (2a + d) (2N-1)$

■ General Ordering Information



INTERDIGITATED MICROSENSOR ELECTRODES (IMES)

IME M - XX50 or 25 - FD, CD or M - *P or U

where M = Au or Pt and XX = 20, 15, 10 or 05 μm

| IME 2050.5 SERIES | GOLD-Au | PLATINUM-Pt | INDIUM TIN OXIDE-ITO |
|-----------------------------------|--------------------|--------------------|----------------------|
| Monolithic un-packaged | IME 2050.5-M-Au | IME 2050.5-M-Pt-U | IME 2050.5-M-ITO-U |
| Monolithic-packaged | IME 2050.5-M-Au-P | IME 2050.5-M-Pt-P | |
| IME 1550.5 SERIES | GOLD-Au | PLATINUM-Pt | INDIUM TIN OXIDE-ITO |
| Monolithic un-packaged | IME 1550.5-M-Au | IME 1550.5-M-Pt-U | IME 1550.5-M-ITO-U |
| Monolithic-packaged | IME 1550.5-M-Au-P | IME 1550.5-M-Pt-P | |
| Combined Differential-un-packaged | IME 1550.5-CD-Au-U | IME 1550.5-CD-Pt-U | IME 1550.5-CD-ITO-U |
| Combined Differential-packaged | IME 1550.5-CD-Au-P | IME 1550.5-CD-Pt-P | |
| Full Differential-un-packaged | IME 1550.5-FD-Au-U | IME 1550.5-FD-Pt-U | IME 1550.5-FD-ITO |
| Full Differential-packaged | IME 1550.5-FD-Au-P | IME 1550.5-FD-Pt-P | |